

INSTRUCTIONS FOR USE

Product Name:

Vesco MultiFil[®] ENFit Oral/Enteral Syringe
Sizes: 1mL, 2.5mL, 3mL, 5mL, 10mL, 20mL and 60mL

Indication for Use:

Vesco MultiFil ENFit Oral/Enteral Syringe is indicated for use as a dispenser, a measuring device, and an oral/enteral fluid transfer device. It is used to deliver fluids into the body orally or enterally. It is intended to be used in clinical or home care settings by users ranging from clinicians to laypersons in all age groups. Vesco MultiFil syringes can be re-used up to six times per day for up to 14 days.

Warnings:

- ENFit Oral/ Enteral syringe has the potential to misconnect with connectors of other healthcare applications.
- The ENFit Oral/Enteral syringe is intended to connect to an Enteral Feeding Tube with an ISO 80369-3 male connector only.
- This syringe can not be used as a storage device.
- The cap (if included) must be twisted to cover all screw threads to ensure a secure fit and/or minimal/no leaking.
- Caps are a choking hazard. Keep out of reach of children.
- Do not modify the syringe. Modification may lead to leakage, inadequate nutrition delivery, or patient harm.
- Do not use if syringe tip is visibly blocked as this may lead to inadequate nutrition delivery.
- The syringe is used for Manual Use Only. Do not use with pumps.

DO NOT REUSE IF:

- Syringe becomes difficult to use. Plunger should easily slide inside barrel.
- Graduation markings become unclear.
- Visible damage to syringe.

Cautions:

- Verify the integrity of each primary container before use.
- Do not use if the package has been previously opened or damaged, or if the device appears damaged.
- Avoid use after the expiration date, for sterile product only.
- **Do not use in a sterile field or environment.**
- Dispose if any signs of damage, or after 14 days.
- Federal law restricts this device to sale by or on the order of a physician.

Instructions:

1. Check the integrity of each primary container before use.
2. Open the package from "PEEL OPEN", and remove the syringe from the package.
3. Remove the syringe's protective cap (if included), and insert the syringe into the liquid. Keep the syringe perpendicular to the liquid level and avoid the syringe touching the walls of the container.
4. Before filling the syringe, the plunger should be in a fully pushed position so that the fluid can be indrawn without air.
5. Draw the liquid by pulling the plunger of the syringe; stop drawing when filling the syringe to the expected volume of the liquid. Ensure the expected graduation line of the scale coincides with the fiducial line of the piston.
6. If you observe air bubbles in the syringe during filling, draw more liquid than needed for the required volume. Keep the syringe tip upward after filling, wait for the air to rise below the nozzle, and slightly push the plunger to expel all air bubbles and excess fluid. Ensure the expected graduation line of the scale coincides with the fiducial line of the piston.
7. For oral feeding, insert the syringe tip into the patient's mouth, then gently push the plunger to deliver the nutrition into the patient's mouth.
8. For use with feeding tube, connect the syringe and feeding tube; push the plunger slowly to deliver the liquid to the feeding tube. Ensure the plunger is pushed to the end of the barrel.

Cleaning & Storage

1. Wash (immediately after each use):
 - Wash hands and ensure clean, clear work surface.
 - Separate syringe and plunger washing both thoroughly in mild soapy water. (Figure B.)
2. Rinse barrel and plunger under cold tap water. (Figure C.)

NOTE: Do not use dishwasher as it may reduce life of syringe.

3. Dry
 - Shake any residual water off syringe components and allow to air dry.
 - Once dry, re-assemble syringe.
4. Store (Figure D.)
 - Store in a clean, dry storage container with a lid.
 - Avoid storage in direct sunlight.
 - Handle with care.
5. Re-Use
 - Use up to six times per day for a maximum of 14 days.

NOTE: Do not use dishcloth or abrasive towel.

Figure B.




Figure C.




Figure A.





Figure D.



Symbol Glossary

Single Patient Multiple Use	Caution Consult Instructions	Not for IV Use	Not made with BPA	Not made with natural rubber latex	Catalogue Number	Non-Sterile
Expiration Date	Electronic Instructions for Use	Keep Dry	Keep Away from Sunlight	Not made with DEHP	Lot Number	Fragile Handle with Care