



Gravity Feeding Sets with ENFit®

- 500 & 1000mL Gravity Bag Sets
- 40mm Screw Cap Closed System Gravity Sets
- Large Bore 1000mL Gravity Bag Set

Features

- Gravity Bag Sets
 - Snap cap with rigid neck for easy filling and closure.
 - Easy-to-read graduations and easy-to-view translucent material.
- Closed System Sets
 - Screw cap allows for flexibility in feeding containers with a 40mm opening.
 - Convenient suspension bag allows for use with smaller containers (baby bottles) for easy hanging.
- Large Bore Gravity Sets
 - Designed to feed thicker formulas faster with 30° tapered bag for quicker, more complete feed.
 - Snap cap with rigid neck for easy filling and closure.
- Roller clamp and 14 drip chamber for easy priming and consistent flow.
- Not made with natural rubber latex, DEHP, Bisphenol A (BPA) or other phthalate plasticizers.

About Vesco Medical

- High Quality Enteral Feeding Devices with ISO 80369-3 compliant ENFit® connectors, ISO 9001 & 13485 Quality System.
- All products made without natural rubber latex, DEHP, Bisphenol A or other phthalate plasticizers
- Competitive Pricing: Savings over competitive brands without compromise in quality
- Unbeatable Service: Friendly, Skilled, & Experienced Customer Care Specialists

Important Safety Information

- To avoid risk of a tubing misconnection, ONLY connect the distal connector to an ENTERAL feeding tube. This device is compatible with an ENFit female connector on the patient's enteral feeding tube.
- Not for IV use.
- This device should be replaced at least every 24 hours.
- Do not store for prolonged periods at extreme temperatures.
- See individual product labeling for directions and important safety information.
- Single-patient use.



Gravity Set Options

Gravity Feed Sets	Item #	Units/Case
40mm Screw Cap Set with ENFit	VED-022	30
500mL Bag Set (LL)* with ENFit	VED-040	30
500mL TOP FILL Bag Set (PVC) with ENFit	VED-042	30
1000mL Bag Set (LL)* with ENFit	VED-034	30
1000mL TOP FILL Bag Set (PVC) with ENFit	VED-046	30
1000mL LARGE BORE Bag Set with ENFit	VED-049	30

*LL Layered Laminate

Setting the Feed Rate

- Vesco Gravity Feeding Sets deliver approximately 14 drops/mL. Adjust roller clamp to control flow rate. Follow chart to set desired rate.

Feeding Rates		
Feeding Rate (mL/Hr)	Drops/Minute	Seconds Between Drops
5	1	52.1
10	3	26
15	4	17.3
20	5	13.9
25	6	10.3
30	7	8.6
50	12	5.1
65	15	4.0
75	18	3.4
100	23	2.6
150	35	1.7
200	47	1.3
250	58	1.0

